




PRICE LIST 2024



Product list 2024

Last revision: February 2024

**All our products are visible and available for purchase on the website
www.ad-maiora.eu**

The images shown in the price list are for illustrative purposes. This means that every single item may not have the endowment shown in the figure below delivery. Variations can be made without prior notice. Ad Maiora S.r.l. is not to be considered responsible for damages deriving from the lack or inaccuracy of the information reported here. The use of the devices described in this catalog is reserved for surgeons who are specifically trained in their application and have completed training courses for the correct implementation of the technique.

Therefore, we invite you to check when ordering with our commercial operators for optimal configurations to meet your specific needs.

For information, quotes and orders:

Via della Costituzione 10, 42025 Cavriago (RE), Italy

Tel. +39 348 868 3311

info@ad-maiora.eu - orders@ad-maiora.eu

Ad Maiora S.r.l.

Ad Maiora aims to develop effective and safe surgical techniques for veterinary patients and create tools that facilitate the surgeon's tasks. Thanks to the intellectual contributions and direct experiences of collaborating surgeons, Ad Maiora sees itself as a constantly evolving entity, striving to provide effective solutions to challenges encountered in surgical practice.

An Ad Maiora tool is a tangible manifestation of the ideas and passion that drive veterinary surgeons. Founded in February 2013 by four individuals, Ad Maiora sought to play a more entrepreneurial role in the design, development, production, and marketing of innovative products for veterinary surgery, leveraging diverse skills.

The fundamental criterion is to address the needs of new techniques, materials, and systems in the rapidly expanding field of veterinary surgery. This involves actively listening to inputs from surgeons and translating them into solutions, thereby serving as the nexus between medical, engineering, and industrial aspects, all while respecting the unique characteristics of the veterinary world.

Collaborations on both the medical and engineering/industrial sides have proliferated, leading to several ongoing development projects, some of which have resulted in patented products.

Index

Ergomed Surgical tables	5
X-Porous TTA	13
TPLO	35
Osti-Lok locked plate system	41
External Fixation	71
Porous titanium cylinders	103
Instruments	105
Tutors	109
K wires and self-tapping threaded pins	113



Ergomed Surgical tables

Ergomed 99

SURGICAL TABLES

Surgical table Ergomed 2M-E



T00010015A

- Up/down and trendelenburg movements
- Single-column structure
- Low voltage control
- Rollover protection device
- Handling wheels
- Can be equipped with all standard accessories
- Top to be added at customer's choice

€ 2,541.00

Surgical table Ergomed 3M-E-18



T00010014A

- Up/down and trendelenburg movements
- Single-column structure
- Lateral tilting
- Low voltage control
- Rollover protection device
- Handling wheels
- Disconnectable cables
- A pair of stabilizers increase stability during orthopedic procedures that will displace large forces off the center of gravity
- Can be equipped with all standard accessories
- Top to be added at customer's choice

€ 3,580.50

Ergomed 99

SURGICAL TABLES

Table top

Stainless steel table top



P00020001A

Stainless steel

€ 142.10

Lexan table top radiolucent



P00020002A

Lexan

€ 179.00

Ergomed table mat



M00080001A

Faux leather

€ 226.40

U mat for Ergomed table



M00080002A

Faux leather

€ 226.40

Ergomed 99

SURGICAL TABLES

Accessories

0.9 meter retention belt



F00010001A

Nylon

€ 19.70

1.8 meter retention belt



F00010002A

Nylon

€ 31.20

Belt hook



G00010003A

Aluminum

€ 13.00

Traction belt with chain



C00030003A

Nylon

€ 24.30

Small dorsal positioner



P00030004A

Stainless steel

€ 53.70

Big dorsal positioner

P00030003A

Stainless steel

€ 80.90

Operating armrest



S000010002A

Stainless steel

€ 92.40

Draping hold-up bar



A00030004A

Stainless steel

€ 98.20

Dripping stand



A00030008A

Stainless steel

€ 63.50

Traction bar



A00030003A

Stainless steel

€ 277.20

Ergomed stabilizing rod



B0001009A

Stainless steel

€ 101.10

Table extension



E00010003A

Lexan

€ 127.10

Universal clamp



M00010015A

Aluminum

€ 133.80

80 mm traction stirrup



S00020001A

Stainless steel

€ 136.90

120 mm traction stirrup



S00020002A

Stainless steel

€ 145.50

160 mm traction stirrup



S00020003A

Stainless steel

€ 159.40

250 mm extension



E0001004A

Stainless steel

€ 86.60

500 mm extension



E0001005A

Stainless steel

€ 92.40

22 kg mechanical dynamometer



S00030006A

Plastic

€ 29.50

Side positioners



S00010001A

Stainless steel

€ 110.00

**Surgical tray
300x400 mm**



S00010004A

Stainless steel

€ 141.00

Accessory rack



C00010002A

Stainless steel

€ 300.30

**Tiltable fluid
collection pan**



V00010006A

Stainless steel

€ 104.00

X-ray cassette tray



S00010003A

Lexan

€ 92.40

Protective padding



P00040001A

Polyurethane

€ 7.00

AD  **MAIORA**®
WE MOVE ON



X-Porous TTA

X-Porous TTA

X-POROUS TTA Titanium wedges 6.3

**Titanium wedge 6.3
3x5**

I00030147A

Titanium

€ 77.00

**Titanium wedge 6.3
3x7**

I00030148A

Titanium

€ 77.00

**Titanium wedge 6.3
4.5x8**



I00030149A

Titanium

€ 77.00

**Titanium wedge 6.3
4.5x11**



I00030150A

Titanium

€ 77.00

**Titanium wedge 6.3
6x11**

I00030151A

Titanium

€ 77.00

**Titanium wedge 6.3
6x14**

I00030152A

Titanium

€ 77.00

**Titanium wedge 6.3
7.5x14**

I00030153A

Titanium

€ 77.00

**Titanium wedge 6.3
7.5x17**

I00030154A

Titanium

€ 77.00

**Titanium wedge 6.3
9x14**



I00030155A

Titanium

€ 77.00

**Titanium wedge 6.3
9x17**



I00030156A

Titanium

€ 77.00

**Titanium wedge 6.3
9x20**



I00030157A

Titanium

€ 77.00

**Titanium wedge 6.3
10.5x17**

I00030158A

Titanium

€ 77.00

**Titanium wedge 6.3
10.5x20**

I00030159A

Titanium

€ 77.00

**Titanium wedge 6.3
12x17**

I00030160A

Titanium

€ 77.00

**Titanium wedge 6.3
12x20**

I00030161A

Titanium

€ 77.00



**Titanium wedge 6.3
12x23**

I00030162A

Titanium

€ 77.00

**Titanium wedge 6.3
13.5x20**

I00030163A

Titanium

€ 77.00



**Titanium wedge 6.3
13.5x23**

I00030164A

Titanium

€ 77.00

**Titanium wedge 6.3
15x21**

I00030165A

Titanium

€ 77.00

**Titanium wedge 6.3
15x24**

I00030166A

Titanium

€ 77.00

X-Porous TTA

X-POROUS TTA

Cortical screws - cross head

**Ø 1.5 x 6 cortical screw,
titanium,
cross head**

V00040051A

Titanium

€10.12

**Ø 1.5 x 8 cortical screw,
titanium,
cross head**

V00040052A

Titanium

€10.12

**Ø 1.5 x 10 cortical screw,
titanium,
cross head**



V00040053A

Titanium

€10.12

**Ø 1.5 x 12 cortical screw,
titanium,
cross head**

V00040054A

Titanium

€10.12

**Ø 1.5 x 14 cortical screw,
titanium,
cross head**

V00040055A

Titanium

€10.12

**Ø 2.0 x 6 cortical screw,
titanium,
cross head**

V00040056A

Titanium

€10.12

**Ø 2.0 x 8 cortical screw,
titanium,
cross head**

V00040040A

Titanium

€ 10.12

**Ø 2.0 x 10 cortical screw,
titanium,
cross head**

V00040041A

Titanium

€ 10.12

**Ø 2.0 x 12 cortical screw,
titanium,
cross head**



V00040042A

Titanium

€ 10.12

**Ø 2.0 x 14 cortical screw,
titanium,
cross head**

V00040057A

Titanium

€ 10.12

**Ø 2.0 x 16 cortical screw,
titanium,
cross head**

V00040058A

Titanium

€ 10.12

X-Porous TTA

X-POROUS TTA

Cortical screws - hexagonal head

**Ø 1.5 x 6 cortical screw,
titanium,
Hex head**

V00040076A

Titanium

€10.58

**Ø 1.5 x 8 cortical screw,
titanium,
Hex head**

V00040077A

Titanium

€10.58

**Ø 1.5 x 10 cortical screw,
titanium,
Hex head**

V00040078A

Titanium

€10.58



**Ø 1.5 x 12 cortical screw,
titanium,
Hex head**

V00040079A

Titanium

€10.58

**Ø 1.5 x 14 cortical screw,
titanium,
Hex head**

V00040080A

Titanium

€10.58

**Ø 2.0 x 6 cortical screw,
titanium,
Hex head**

V00040081A

Titanium

€10.58

**Ø 2.0 x 8 cortical screw,
titanium,
Hex head**

V00040082A

Titanium

€ 10.58

**Ø 2.0 x 10 cortical screw,
titanium,
Hex head**

V00040083A

Titanium

€ 10.58

**Ø 2.0 x 12 cortical screw
titanium,
Hex head**



V00040084A

Titanium

€ 10.58

**Ø 2.0 x 14 cortical screw,
titanium,
Hex head**

V00040085A

Titanium

€ 10.58

**Ø 2.0 x 16 cortical screw,
titanium,
Hex head**

V00040086A

Titanium

€ 10.58

X-Porous TTA

X-POROUS TTA

Cortical screws - hexagonal head

**Ø 2.4 x 10 cortical screw,
titanium,
Hex head**

V00040043A

Titanium

€ 10.12

**Ø 2.4 x 12 cortical screw,
titanium,
Hex head**

V00040039A

Titanium

€ 10.12

**Ø 2.4 x 14 cortical screw,
titanium,
Hex head**

V00040044A

Titanium

€ 10.12



**Ø 2.4 x 16 cortical screw,
titanium,
Hex head**

V00040059A

Titanium

€ 10.12

**Ø 2.4 x 18 cortical screw,
titanium,
Hex head**

V00040060A

Titanium

€ 10.12

**Ø 2.4 x 20 cortical screw,
titanium,
Hex head**

V00040061A

Titanium

€ 10.12

**Ø 2.4 x 22 cortical screw,
titanium,
Hex head**

V00040062A

Titanium

€ 10.12

**Ø 2.4 x 24 cortical screw,
titanium,
Hex head**

V00040063A

Titanium

€ 10.12

**Ø 2.4 x 26 cortical screw,
titanium,
Hex head**



V00040064A

Titanium

€ 10.12

**Ø 2.4 x 28 cortical screw,
titanium,
Hex head**

V00040065A

Titanium

€ 10.12

**Ø 2.7 x 14 cortical screw,
titanium,
Hex head**

V00040066A

Titanium

€ 10.12

**Ø 2.7 x 16 cortical screw,
titanium,
Hex head**



V00040067A

Titanium

€ 10.12

**Ø 2.7 x 18 cortical screw,
titanium,
Hex head**

V00040068A

Titanium

€ 10.12

**Ø 2.7 x 20 cortical screw,
titanium,
Hex head**

V00040045A

Titanium

€ 10.12

**Ø 2.7 x 22 cortical screw,
titanium,
Hex head**

V00040069A

Titanium

€ 10.12

**Ø 2.7 x 24 cortical screw,
titanium,
Hex head**



V00040046A

Titanium

€ 10.12

**Ø 2.7 x 26 cortical screw,
titanium,
Hex head**

V00040047A

Titanium

€ 10.12

**Ø 2.7 x 28 cortical screw,
titanium,
Hex head**

V00040075A

Titanium

€ 10.12

**Ø 3.5 x 16 cortical screw,
titanium,
Hex head**

V00040048A

Titanium

€ 10.12

**Ø 3.5 x 18 cortical screw,
titanium,
Hex head**



V00040049A

Titanium

€ 10.12

**Ø 3.5 x 20 cortical screw,
titanium,
Hex head**

V00040070A

Titanium

€ 10.12

**Ø 3.5 x 22 cortical screw,
titanium,
Hex head**

V00040050A

Titanium

€ 10.12

**Ø 3.5 x 24 cortical screw,
titanium,
Hex head**

V00040071A

Titanium

€ 10.12

**Ø 3.5 x 26 cortical screw,
titanium,
Hex head**



V00040072A

Titanium

€ 10.12

**Ø 3.5 x 28 cortical screw,
titanium,
Hex head**

V00040073A

Titanium

€ 10.12

**Ø 3.5 x 30 cortical screw,
titanium,
Hex head**

V00040074A

Titanium

€ 10.12

X-Porous TTA

X-POROUS TTA

Locking screws - Hexagonal head

**Ø 2.0 x 10 locking screw,
titanium,
Hex head**

V00050010A

Titanium

€17.60

**Ø 2.0 x 12 locking screw,
titanium,
Hex head**

V00050009A

Titanium

€17.60

**Ø 2.0 x 14 locking screw,
titanium,
Hex head**



V00050011A

Titanium

€17.60

**Ø 2.0 x 16 locking screw,
titanium,
Hex head**

V00050012A

Titanium

€17.60

**Ø 2.0 x 18 locking screw,
titanium,
Hex head**

V00050095A

Titanium

€17.60

**Ø 2.0 x 20 locking screw,
titanium,
Hex head**

V00050096A

Titanium

€17.60

**Ø 2.4 x 10 locking screw,
titanium,
Hex head**

V00050007A

Titanium

€ 17.60

**Ø 2.4 x 12 locking screw,
titanium,
Hex head**

V00050001A

Titanium

€ 17.60

**Ø 2.4 x 14 locking screw,
titanium,
Hex head**

V00050002A

Titanium

€ 17.60

**Ø 2.4 x 16 locking screw,
titanium,
Hex head**

V00050003A

Titanium

€ 17.60

**Ø 2.4 x 18 locking screw,
titanium,
Hex head**

V00050097A

Titanium

€ 17.60

**Ø 2.4 x 20 locking screw,
titanium,
Hex head**

V00050098A

Titanium

€ 17.60

**Ø 2.7 x 10 locking screw,
titanium,
Hex head**

V00050008A

Titanium

€ 17.60



**Ø 2.7 x 12 locking screw,
titanium,
Hex head**

V00050004A

Titanium

€ 17.60

**Ø 2.7 x 14 locking screw,
titanium,
Hex head**

V00050005A

Titanium

€ 17.60

**Ø 2.7 x 16 locking screw,
titanium,
Hex head**



V00050006A

Titanium

€ 17.60

**Ø 2.7 x 18 locking screw,
titanium,
Hex head**

V00050099A

Titanium

€ 17.60

**Ø 2.7 x 20 locking screw,
titanium,
Hex head**

V00050100A

Titanium

€ 17.60

X-Porous TTA

X-POROUS TTA Titanium locked plates

**Titanium locked plate
7S3.0B Right**



I00020052A

Titanium

€ 38.50

**Titanium locked plate
7S3.0B Left**



I00020057A

Titanium

€ 38.50

**Titanium locked plate
7M.0B Right**



I00020051A

Titanium

€ 38.50

**Titanium locked plate
7M3.0B Left**



I00020056A

Titanium

€ 38.50

**Titanium locked plate
7L3.0B Right**



I00020050A

Titanium

€ 38.50

**Titanium locked plate
7L3.0B Left**

I00020055A

Titanium

€ 38.50

**Titanium locked plate
7SL3.0B Right**

I00020075A

Titanium

€ 38.50

**Titanium locked plate
7SL3.0B Left**



I00020076A

Titanium

€ 38.50

**Titanium locked plate
83S.0B Right**

I00020054A

Titanium

€ 38.50

**Titanium locked plate
83S.0B Left**



I00020059A

Titanium

€ 38.50

**Titanium locked plate
8L3.0B Right**

I00020053A

Titanium

€ 38.50

**Titanium locked plate
8L3.0B Left**



I00020058A

Titanium

€ 38.50

**Titanium locked plate
8SL3.0B Right**

I00020091A

Titanium

€ 38.50

**Titanium locked plate
8SL3.0B Left**



I00020092A

Titanium

€ 38.50

**Titanium locked plate
8SM3.0B Right**



I00020093A

Titanium

€ 38.50

**Titanium locked plate
8SM3.0B Left**



I00020094A

Titanium

€ 38.50

X-Porous TTA

X-POROUS TTA Titanium plates

Titanium plate 4S3.0



I00020060E

Titanium

€ 25.30

Titanium plate 4M3.0



I00020061E

Titanium

€ 25.30

Titanium plate 4L3.0



I00020062E

Titanium

€ 25.30

Titanium plate 7S3.0



I00020047E

Titanium

€ 25.30

Titanium plate 7M3.0



I00020046E

Titanium

€ 25.30

Titanium plate 7L3.0



I00020045E

Titanium

€ 25.30

Titanium plate 7SL3.0



I00020077E

Titanium

€ 25.30

Titanium plate 8S3.0



I00020049E

Titanium

€ 25.30

Titanium plate 8SM3.0



I00020089A

Titanium

€ 25.30

Titanium plate 8SL3.0



I00020090A

Titanium

€ 25.30

Titanium plate 8L3.0



I00020048E

Titanium

€ 25.30

X-Porous TTA

X-POROUS TTA Containers

**Container for implants
and instruments for
X-Porous TTA**



C00040021A

Thermoplastic material

€ 220.00

**Extended container for
X-Porous TTA implants**



C00040027A

Technopolymer

€ 330.00

AD  **MAIORA**®
WE MOVE ON



TPLO

Locked plate system

TPLO

Standard TPLO Set

Standard TPLO Set

This configuration is ideal for sites focused on TPLO procedures only.

€ 5,651.94

Kit composition

Code	Product description	Quantity
I00020095A	TPLO Plate Large, Right	1
I00020096A	TPLO Plate Large, Left	1
I00020106A	TPLO Plate Large, Flat	1
I00020107A	TPLO Plate Medium, Right	1
I00020108A	TPLO Plate Medium, Left	1
I00020109A	TPLO Plate Medium, Flat	1
I00020110A	TPLO Plate Small, Right	1
I00020111A	TPLO Plate Small, Left	1
I00020112A	TPLO Plate Small, Flat	2
V00050075A	Screw Ø 3.5 x 10 locked, titanium, torx head	4
V00050076A	Screw Ø 3.5 x 12 locked, titanium, torx head	4
V00050077A	Screw Ø 3.5 x 14 locked, titanium, torx head	6
V00050078A	Screw Ø 3.5 x 16 locked, titanium, torx head	6
V00050079A	Screw Ø 3.5 x 18 locked, titanium, torx head	6
V00050080A	Screw Ø 3.5 x 20 locked, titanium, torx head	6
V00050081A	Screw Ø 3.5 x 22 locked, titanium, torx head	4
V00050082A	Screw Ø 3.5 x 24 locked, titanium, torx head	4
V00050083A	Screw Ø 3.5 x 26 locked, titanium, torx head	2
V00050084A	Screw Ø 3.5 x 28 locked, titanium, torx head	2
V00050085A	Screw Ø 3.5 x 30 locked, titanium, torx head	2
V00050086A	Screw Ø 3.5 x 32 locked, titanium, torx head	2
V00050087A	Screw Ø 3.5 x 34 locked, titanium, torx head	2
V00050088A	Screw Ø 3.5 x 36 locked, titanium, torx head	2
V00050089A	Screw Ø 3.5 x 38 locked, titanium, torx head	2
V00050090A	Screw Ø 3.5 x 40 locked, titanium, torx head	4
V00050064A	Screw Ø 2.7 x 8 locked, titanium, torx head	4
V00050065A	Screw Ø 2.7 x 10 locked, titanium, torx head	6
V00050066A	Screw Ø 2.7 x 12 locked, titanium, torx head	6
V00050067A	Screw Ø 2.7 x 14 locked, titanium, torx head	6
V00050068A	Screw Ø 2.7 x 16 locked, titanium, torx head	4
V00050069A	Screw Ø 2.7 x 18 locked, titanium, torx head	2
V00050070A	Screw Ø 2.7 x 20 locked, titanium, torx head	2
V00050071A	Screw Ø 2.7 x 22 locked, titanium, torx head	2
V00050072A	Screw Ø 2.7 x 24 locked, titanium, torx head	2
V00040135A	Screw Ø 3.5 x 10 cortical, titanium, torx head	2
V00040136A	Screw Ø 3.5 x 12 cortical, titanium, torx head	2
V00040137A	Screw Ø 3.5 x 14 cortical, titanium, torx head	2
V00040138A	Screw Ø 3.5 x 16 cortical, titanium, torx head	2
V00040139A	Screw Ø 3.5 x 18 cortical, titanium, torx head	2
V00040140A	Screw Ø 3.5 x 20 cortical, titanium, torx head	2
V00040141A	Screw Ø 3.5 x 22 cortical, titanium, torx head	2
V00040142A	Screw Ø 3.5 x 24 cortical, titanium, torx head	2

TPLO

Standard TPLO Set

Standard TPLO Set

This configuration is ideal for sites focused on TPLO procedures only.

€ 5,651.94

Kit composition

Code	Product description	Quantity
V00040143A	Screw Ø 3.5 x 26 cortical, titanium, torx head	2
V00040144A	Screw Ø 3.5 x 28 cortical, titanium, torx head	2
V00040145A	Screw Ø 3.5 x 30 cortical, titanium, torx head	2
V00040113A	Screw Ø 2.7 x 10 cortical, titanium, torx head	2
V00040114A	Screw Ø 2.7 x 12 cortical, titanium, torx head	2
V00040115A	Screw Ø 2.7 x 14 cortical, titanium, torx head	2
V00040116A	Screw Ø 2.7 x 16 cortical, titanium, torx head	2
V00040117A	Screw Ø 2.7 x 18 cortical, titanium, torx head	2
V00040118A	Screw Ø 2.7 x 20 cortical, titanium, torx head	2
V00040119A	Screw Ø 2.7 x 22 cortical, titanium, torx head	2
V00040120A	Screw Ø 2.7 x 24 cortical, titanium, torx head	2
V00040121A	Screw Ø 2.7 x 26 cortical, titanium, torx head	2
V00040122A	Screw Ø 2.7 x 28 cortical, titanium, torx head	2
V00040123A	Screw Ø 2.7 x 30 cortical, titanium, torx head	2
V00040124A	Screw Ø 2.7 x 32 cortical, titanium, torx head	2
V00040125A	Screw Ø 2.7 x 34 cortical, titanium, torx head	2
V00040126A	Screw Ø 2.7 x 36 cortical, titanium, torx head	2
V00040127A	Screw Ø 2.7 x 38 cortical, titanium, torx head	2
V00040128A	Screw Ø 2.7 x 40 cortical, titanium, torx head	2
P00060006A	3,5-mm drill bit with quick coupling	1
P00060007A	2.7-mm drill bit with quick coupling	1
U00010038A	Screwdriver T-10	1
S00050003A	Drill guide with lock	2
P00090003A	In-plane bender	2
U00010042A	Handle grip nut	4
P00040009A	Angled grip nut, small	1
P00040010A	Angled grip nut, medium	1
P00040012A	Angled grip nut, large	1
P00040006A	Grip nut, small	2
P00040007A	Grip nut, medium	2
P00040008A	Grip nut, large	2
M00070002A	Depth gauge, 40 mm	1
U00010039A	Bender direct connection	1
U00010040A	AO driver insert, T10	1
C00040029A	Osti-lok case	1

Locked plate system

TPLO

TPLO Plates

**TPLO Plate Large,
Right**

I00020095A

Titanium

€ 123.20

**TPLO Plate Large,
Left**



I00020096A

Titanium

€ 123.20

**TPLO Plate Large,
Flat**

I00020106A

Titanium

€ 123.20

**TPLO Plate Medium,
Right**

I00020107A

Titanium

€ 118.80

**TPLO Plate Medium,
Left**



I00020108A

Titanium

€ 118.80

**TPLO Plate Medium,
Flat**

I00020109A

Titanium

€ 118.80

**PLO Plate Small,
Right**

I00020110A

Titanium

€ 104.50

**PLO Plate Small,
Left**



I00020111A

Titanium

€ 104.50

**PLO Plate Small,
Flat**

I00020112A

Titanium

€ 104.50

AD  **MAIORA**®
WE MOVE ON



Osti-Lok locked plate system

Osti-Lok locked plate system

OSTI-LOK

Standard Osti-Lok Set

Standard Osti-Lok Set



Using the latest technologies, Osti-Lok was conceived from the ground up for maximal implant versatility, while requiring minimal instrumentation.

€ 7,014.04

Kit composition

Code	Product description	Quantity
I00020063A	Straight titanium small plate PS-100	1
I00020064A	Straight titanium medium plate PM-100	1
I00020065A	Straight titanium large plate PL-100	1
I00020066A	Titanium small cross plate CPS-100	1
I00020067A	Titanium medium cross plate CPM-100	1
I00020068A	Titanium large cross plate CPL-100	1
I00020078A	Small plate, high strength PSHS-100	1
I00020079A	Medium plate, high strength PMHS-100	1
I00020080A	Large plate, high strength PLHS-100	1
I00020081A	Small plate, 16 hole PS16-100	1
I00020082A	Medium plate, 16 hole PM16-100	1
I00020083A	Large plate, 16 hole PL16-100	1
V00050075A	Screw Ø 3.5 x 10 locked, titanium, torx head	2
V00050076A	Screw Ø 3.5 x 12 locked, titanium, torx head	4
V00050077A	Screw Ø 3.5 x 14 locked, titanium, torx head	4
V00050078A	Screw Ø 3.5 x 16 locked, titanium, torx head	6
V00050079A	Screw Ø 3.5 x 18 locked, titanium, torx head	6
V00050080A	Screw Ø 3.5 x 20 locked, titanium, torx head	6
V00050081A	Screw Ø 3.5 x 22 locked, titanium, torx head	6
V00050082A	Screw Ø 3.5 x 24 locked, titanium, torx head	4
V00050083A	Screw Ø 3.5 x 26 locked, titanium, torx head	4
V00050084A	Screw Ø 3.5 x 28 locked, titanium, torx head	2
V00050085A	Screw Ø 3.5 x 30 locked, titanium, torx head	2
V00050086A	Screw Ø 3.5 x 32 locked, titanium, torx head	2
V00050087A	Screw Ø 3.5 x 34 locked, titanium, torx head	2
V00050088A	Screw Ø 3.5 x 36 locked, titanium, torx head	2
V00050089A	Screw Ø 3.5 x 38 locked, titanium, torx head	2
V00050090A	Screw Ø 3.5 x 40 locked, titanium, torx head	2
V00050064A	Screw Ø 2.7 x 8 locked, titanium, torx head	4
V00050065A	Screw Ø 2.7 x 10 locked, titanium, torx head	4
V00050066A	Screw Ø 2.7 x 12 locked, titanium, torx head	6
V00050067A	Screw Ø 2.7 x 14 locked, titanium, torx head	6
V00050068A	Screw Ø 2.7 x 16 locked, titanium, torx head	6
V00050069A	Screw Ø 2.7 x 18 locked, titanium, torx head	4
V00050070A	Screw Ø 2.7 x 20 locked, titanium, torx head	2
V00050071A	Screw Ø 2.7 x 22 locked, titanium, torx head	2
V00050072A	Screw Ø 2.7 x 24 locked, titanium, torx head	2
V00040135A	Screw Ø 3.5 x 10 cortical, titanium, torx head	2
V00040136A	Screw Ø 3.5 x 12 cortical, titanium, torx head	2
V00040137A	Screw Ø 3.5 x 14 cortical, titanium, torx head	2
V00040138A	Screw Ø 3.5 x 16 cortical, titanium, torx head	2
V00040139A	Screw Ø 3.5 x 18 cortical, titanium, torx head	2
V00040140A	Screw Ø 3.5 x 20 cortical, titanium, torx head	2

Osti-Lok locked plate system

OSTI-LOK

Standard Osti-Lok Set

Standard Osti-Lok Set



Using the latest technologies, Osti-Lok was conceived from the ground up for maximal implant versatility, while requiring minimal instrumentation.

€ 7,014.04

Kit composition

Code	Product description	Quantity
V00040141A	Screw Ø 3.5 x 22 cortical, titanium, torx head	2
V00040142A	Screw Ø 3.5 x 24 cortical, titanium, torx head	2
V00040143A	Screw Ø 3.5 x 26 cortical, titanium, torx head	2
V00040144A	Screw Ø 3.5 x 28 cortical, titanium, torx head	2
V00040145A	Screw Ø 3.5 x 30 cortical, titanium, torx head	2
V00040113A	Screw Ø 2.7 x 10 cortical, titanium, torx head	2
V00040114A	Screw Ø 2.7 x 12 cortical, titanium, torx head	2
V00040115A	Screw Ø 2.7 x 14 cortical, titanium, torx head	2
V00040116A	Screw Ø 2.7 x 16 cortical, titanium, torx head	2
V00040117A	Screw Ø 2.7 x 18 cortical, titanium, torx head	2
V00040118A	Screw Ø 2.7 x 20 cortical, titanium, torx head	2
V00040119A	Screw Ø 2.7 x 22 cortical, titanium, torx head	2
V00040120A	Screw Ø 2.7 x 24 cortical, titanium, torx head	2
V00040121A	Screw Ø 2.7 x 26 cortical, titanium, torx head	2
V00040122A	Screw Ø 2.7 x 28 cortical, titanium, torx head	2
V00040123A	Screw Ø 2.7 x 30 cortical, titanium, torx head	2
V00040124A	Screw Ø 2.7 x 32 cortical, titanium, torx head	2
V00040125A	Screw Ø 2.7 x 34 cortical, titanium, torx head	2
V00040126A	Screw Ø 2.7 x 36 cortical, titanium, torx head	2
V00040127A	Screw Ø 2.7 x 38 cortical, titanium, torx head	2
V00040128A	Screw Ø 2.7 x 40 cortical, titanium, torx head	2
P00060006A	3,5-mm drill bit with quick coupling	1
P00060007A	2.7-mm drill bit with quick coupling	1
U00010038A	Screwdriver T-10	1
S00050003A	Drill guide with lock	2
P00090003A	In-plane bender	2
U00010042A	Handle grip nut	4
P00040009A	Angled grip nut, small	1
P00040010A	Angled grip nut, medium	1
P00040012A	Angled grip nut, large	1
P00040006A	Grip nut, small	2
P00040007A	Grip nut, medium	2
P00040008A	Grip nut, large	2
M00070002A	Depth gauge, 40 mm	1
U00010039A	Bender direct connection	1
U00010040A	AO driver insert, T10	1
C00040029A	Osti-lok case	1

Osti-Lok locked plate system

OSTI-LOK

Basic Osti-Lok Set

Basic Osti-Lok Set



Using the latest technologies, Osti-Lok was conceived from the ground up for maximal implant versatility, while requiring minimal instrumentation.

€ 5,981.58

Kit composition

Code	Product description	Quantity
I00020063A	Straight titanium small plate PS-100	1
I00020064A	Straight titanium medium plate PM-100	1
I00020065A	Straight titanium large plate PL-100	1
I00020066A	Titanium small cross plate CPS-100	1
I00020067A	Titanium medium cross plate CPM-100	1
I00020068A	Titanium large cross plate CPL-100	1
I00020078A	Small plate, high strength PSHS-100	1
I00020079A	Medium plate, high strength PMHS-100	1
I00020080A	Large plate, high strength PLHS-100	1
I00020081A	Small plate, 16 hole PS16-100	1
I00020082A	Medium plate, 16 hole PM16-100	1
I00020083A	Large plate, 16 hole PL16-100	1
V00050075A	Screw Ø 3.5 x 10 locked, titanium, torx head	2
V00050076A	Screw Ø 3.5 x 12 locked, titanium, torx head	4
V00050077A	Screw Ø 3.5 x 14 locked, titanium, torx head	4
V00050078A	Screw Ø 3.5 x 16 locked, titanium, torx head	6
V00050079A	Screw Ø 3.5 x 18 locked, titanium, torx head	6
V00050080A	Screw Ø 3.5 x 20 locked, titanium, torx head	6
V00050081A	Screw Ø 3.5 x 22 locked, titanium, torx head	6
V00050082A	Screw Ø 3.5 x 24 locked, titanium, torx head	4
V00050083A	Screw Ø 3.5 x 26 locked, titanium, torx head	4
V00050084A	Screw Ø 3.5 x 28 locked, titanium, torx head	2
V00050085A	Screw Ø 3.5 x 30 locked, titanium, torx head	2
V00050086A	Screw Ø 3.5 x 32 locked, titanium, torx head	2
V00050087A	Screw Ø 3.5 x 34 locked, titanium, torx head	2
V00050088A	Screw Ø 3.5 x 36 locked, titanium, torx head	2
V00050089A	Screw Ø 3.5 x 38 locked, titanium, torx head	2
V00050090A	Screw Ø 3.5 x 40 locked, titanium, torx head	2
V00050064A	Screw Ø 2.7 x 8 locked, titanium, torx head	4
V00050065A	Screw Ø 2.7 x 10 locked, titanium, torx head	4
V00050066A	Screw Ø 2.7 x 12 locked, titanium, torx head	6
V00050067A	Screw Ø 2.7 x 14 locked, titanium, torx head	6
V00050068A	Screw Ø 2.7 x 16 locked, titanium, torx head	6
V00050069A	Screw Ø 2.7 x 18 locked, titanium, torx head	4
V00050070A	Screw Ø 2.7 x 20 locked, titanium, torx head	2
V00050071A	Screw Ø 2.7 x 22 locked, titanium, torx head	2
V00050072A	Screw Ø 2.7 x 24 locked, titanium, torx head	2
V00040066A	Screw Ø 2.7 x 14 cortical, titanium, hex head	2
V00040067A	Screw Ø 2.7 x 16 cortical, titanium, hex head	2
V00040068A	Screw Ø 2.7 x 18 cortical, titanium, hex head	2
V00040045A	Screw Ø 2.7 x 20 cortical, titanium, hex head	2
V00040069A	Screw Ø 2.7 x 22 cortical, titanium, hex head	2

Osti-Lok locked plate system

OSTI-LOK

Basic Osti-Lok Set

Basic Osti-Lok Set



Using the latest technologies, Osti-Lok was conceived from the ground up for maximal implant versatility, while requiring minimal instrumentation.

€ 5,981.58

Kit composition

Code	Product description	Quantity
V00040046A	Screw Ø 2.7 x 24 cortical, titanium, hex head	2
V00040047A	Screw Ø 2.7 x 26 cortical, titanium, hex head	2
V00040075A	Screw Ø 2.7 x 28 cortical, titanium, hex head	2
P00060006A	3,5-mm drill bit with quick coupling	2
P00060007A	2.7-mm drill bit with quick coupling	2
U00010038A	Screwdriver T-10	1
U00010056A	Fiber screwdriver for hexagonal screws Ø 2.5 mm	1
S00050003A	Drill guide with lock	2
P00090003A	In-plane bender	2
U00010042A	Handle grip nut	4
P00040009A	Angled grip nut, small	1
P00040010A	Angled grip nut, medium	1
P00040012A	Angled grip nut, large	1
P00040006A	Grip nut, small	2
P00040007A	Grip nut, medium	2
P00040008A	Grip nut, large	2
M00070002A	Depth gauge, 40 mm	1
U00010039A	Bender direct connection	1
U00010040A	AO driver insert, T10	1
C00040029A	Osti-lok case	1

Osti-Lok locked plate system

OSTI-LOK

Locking plates

**Straight small plate,
13 holes**

I00020063A

Titanium

€ 145.20

**Straight medium plate,
13 holes**



I00020064A

Titanium

€ 175.70

**Straight large plate,
13 holes**

I00020065A

Titanium

€ 202.40

**Small cross plate,
13 holes**

I00020066A

Titanium

€ 183.70

**Medium cross plate,
13 holes**



I00020067A

Titanium

€ 204.70

**Large cross plate,
13 holes**

I00020068A

Titanium

€ 217.80

**Small high strength plate,
10 holes**

I00020078A

Titanium

€ 167.20

**Medium high strength plate,
10 holes**

I00020079A

Titanium

€ 189.20

**Large high strength plate,
10 holes**

I00020080A

Titanium

€ 228.80

**Straight small plate,
16 hole**

I00020081A

Titanium

€ 184.80

**Straight medium plate,
16 hole**

I00020082A

Titanium

€ 215.60

**Straight large plate,
16 hole**

I00020083A

Titanium

€ 248.60



Osti-Lok locked plate system

OSTI-LOK

Osti-Lok locking screws

**Ø 3.5 x 10 locking screw,
titanium,
torx head**

V00050075A

Titanium

€ 23.76

**Ø 3.5 x 12 locking screw,
titanium,
torx head**

V00050076A

Titanium

€ 23.76

**Ø 3.5 x 14 locking screw,
titanium,
torx head**



V00050077A

Titanium

€ 23.76

**Ø 3.5 x 16 locking screw,
titanium,
torx head**

V00050078A

Titanium

€ 23.76

**Ø 3.5 x 18 locking screw,
titanium,
torx head**

V00050079A

Titanium

€ 23.76

**Ø 3.5 x 20 locking screw,
titanium,
torx head**

V00050080A

Titanium

€ 23.76

**Ø 3.5 x 22 locking screw,
titanium,
torx head**

V00050081A

Titanium

€ 23.76

**Ø 3.5 x 24 locking screw,
titanium,
torx head**

V00050082A

Titanium

€ 23.76

**Ø 3.5 x 26 locking screw,
titanium,
torx head**

V00050083A

Titanium

€ 23.76

**Ø 3.5 x 28 locking screw,
titanium,
torx head**



V00050084A

Titanium

€ 23.76

**Ø 3.5 x 30 locking screw,
titanium,
torx head**

V00050085A

Titanium

€ 23.76

**Ø 3.5 x 32 locking screw,
titanium,
torx head**

V00050086A

Titanium

€ 23.76

**Ø 3.5 x 34 locking screw,
titanium,
torx head**

V00050087A

Titanio

€ 23.76

**Ø 3.5 x 36 locking screw,
titanium,
torx head**

V00050088A

Titanium

€ 23.76

**Ø 3.5 x 38 locking screw,
titanium,
torx head**



V00050089A

Titanium

€ 23.76

**Ø 3.5 x 40 locking screw,
titanium,
torx head**

V00050090A

Titanium

€ 23.76

**Ø 2.7 x 8 locking screw,
titanium,
torx head**

V00050064A

Titanium

€ 23.76

**Ø 2.7 x 10 locking screw,
titanium,
torx head**



V00050065A

Titanium

€ 23.76

**Ø 2.7 x 12 locking screw,
titanium,
torx head**

V00050066A

Titanium

€ 23.76

**Ø 2.7 x 14 locking screw,
titanium,
torx head**

V00050067A

Titanium

€ 23.76

**Ø 2.7 x 16 locking screw,
titanium,
torx head**

V00050068A

Titanium

€ 23.76

**Ø 2.7 x 18 locking screw,
titanium,
torx head**

V00050069A

Titanium

€ 23.76

**Ø 2.7 x 20 locking screw,
titanium,
torx head**



V00050070A

Titanium

€ 23.76

**Ø 2.7 x 22 locking screw,
titanium,
torx head**

V00050071A

Titanium

€ 23.76

**Ø 2.7 x 24 locking screw,
titanium,
torx head**

V00050072A

Titanium

€ 23.76

Osti-Lok locked plate system

OSTI-LOK

Osti-Lok cortical screws

**Ø 3.5 x 10 cortical screw,
titanium,
torx head**

V00040135A

Titanium

€ 20.14

**Ø 3.5 x 12 cortical screw,
titanium,
torx head**

V00040136A

Titanium

€ 20.14

**Ø 3.5 x 14 cortical screw,
titanium,
torx head**



V00040137A

Titanium

€ 20.14

**Ø 3.5 x 16 cortical screw,
titanium,
torx head**

V00040138A

Titanium

€ 20.14

**Ø 3.5 x 18 cortical screw,
titanium,
torx head**

V00040139A

Titanium

€ 20.14

**Ø 3.5 x 20 cortical screw,
titanium,
torx head**

V00040140A

Titanium

€ 20.14

**Ø 3.5 x 22 cortical screw,
titanium,
torx head**

V00040141A

Titanium

€ 20.14

**Ø 3.5 x 24 cortical screw,
titanium,
torx head**

V00040142A

Titanium

€ 20.14

**Ø 3.5 x 26 cortical screw,
titanium,
torx head**



V00040143A

Titanium

€ 20.14

**Ø 3.5 x 28 cortical screw,
titanium,
torx head**

V00040144A

Titanium

€ 20.14

**Ø 3.5 x 30 cortical screw,
titanium,
torx head**

V00040145A

Titanium

€ 20.14

**Ø 2.7 x 10 cortical screw,
titanium,
torx head**

V00040113A

Titanium

€ 20.14

**Ø 2.7 x 12 cortical screw,
titanium,
torx head**



V00040114A

Titanium

€ 20.14

**Ø 2.7 x 14 cortical screw,
titanium,
torx head**

V00040115A

Titanium

€ 20.14

**Ø 2.7 x 16 cortical screw,
titanium,
torx head**

V00040116A

Titanium

€ 20.14

**Ø 2.7 x 18 cortical screw,
titanium,
torx head**

V00040117A

Titanium

€ 20.14

**Ø 2.7 x 20 cortical screw,
titanium,
torx head**



V00040118A

Titanium

€ 20.14

**Ø 2.7 x 22 cortical screw,
titanium,
torx head**

V00040119A

Titanium

€ 20.14

**Ø 2.7 x 24 cortical screw,
titanium,
torx head**

V00040120A

Titanium

€ 20.14

**Ø 2.7 x 26 cortical screw,
titanium,
torx head**

V00040121A

Titanium

€ 20.14

**Ø 2.7 x 28 cortical screw,
titanium,
torx head**

V00040122A

Titanium

€ 20.14

**Ø 2.7 x 30 cortical screw,
titanium,
torx head**

V00040123A

Titanium

€ 20.14

**Ø 2.7 x 32 cortical screw,
titanium,
torx head**

V00040124A

Titanium

€ 20.14

**Ø 2.7 x 34 cortical screw,
titanium,
torx head**



V00040125A

Titanium

€ 20.14

**Ø 2.7 x 36 cortical screw,
titanium,
torx head**

V00040126A

Titanium

€ 20.14

**Ø 2.7 x 38 cortical screw,
titanium,
torx head**

V00040127A

Titanium

€ 20.14

**Ø 2.7 x 40 cortical screw,
titanium,
torx head**

V00040128A

Titanium

€ 20.14

Osti-Lok locked plate system

OSTI-LOK

Osti-Lok replacement elements

**Angled grip nut,
small**



P00040009A

Titanium

€ 44.00

**Angled grip nut,
medium**



P00040010A

Titanium

€ 44.00

**Angled grip nut,
large**



P00040012A

Titanium

€ 44.00

**Grip nut,
small**



P00040006A

Titanium

€ 44.00

**Grip nut,
medium**



P00040007A

Titanium

€ 44.00

**Grip nut,
large**



P00040008A

Titanium

€ 44.00

Ø 3.5-mm drill bit with quick coupling



P00060006A

Titanium

€ 66.00

Ø 2.7 mm drill bit with quick coupling



P00060007A

Titanium

€ 66.00

Screwdriver T10



U00010038A

Titanium

€ 88.00

AO driver insert, T10



U00010040A

Titanium

€ 44.00

Drill guide with lock



S0005007A

Titanium

€ 44.00

In-plane bender



P00090003A

Titanium

€ 44.00

Depth gauge, 40 mm



M00070002A

Titanium

€ 77.00

Handle grip nut



U00010042A

Titanium

€ 44.00

Bender direct connection



U00010039A

Titanium

€ 44.00

Osti-Lok case



C00040029A

Titanium

€ 451.00

Osti-Lok locked plate system

OSTI-LOK MINI

Osti-Lok mini set

Osti-Lok mini set



Suited for the treatment of patients up to 10 kg

€ 5,270.54

Kit composition

Code	Product description	Quantity
I00020069A	Cross plate CPX 20-050	1
I00020071A	Straight plates PXL 20-050	1
I00020084A	High strength plate, PX20HS-050	1
I00020085A	Radius plate Y - PRX20-050	1
I00020072A	Cross plate CPX 15-050	1
I00020074A	Straight plate PXL 15-050	1
I00020086A	High strength plate, PX15HS-050	1
I00020087A	Radius plate, PRX15-050	1
V00050073A	Screw Ø 2.0 x 6 locked, titanium, square head	4
V00050038A	Screw Ø 2.0 x 8 locked, titanium, square head	4
V00050039A	Screw Ø 2.0 x 10 locked, titanium, square head	4
V00050040A	Screw Ø 2.0 x 12 locked, titanium, square head	4
V00050041A	Screw Ø 2.0 x 14 locked, titanium, square head	4
V00050042A	Screw Ø 2.0 x 16 locked, titanium, square head	4
V00050093A	Screw Ø 2.0 x 18 locked, titanium, square head	4
V00050094A	Screw Ø 2.0 x 20 locked, titanium, square head	4
V00050074A	Screw Ø 1.5 x 6 locked, titanium, square head	4
V00050043A	Screw Ø 1.5 x 8 locked, titanium, square head	4
V00050044A	Screw Ø 1.5 x 10 locked, titanium, square head	4
V00050045A	Screw Ø 1.5 x 12 locked, titanium, square head	4
V00050046A	Screw Ø 1.5 x 14 locked, titanium, square head	4
V00050047A	Screw Ø 1.5 x 16 locked, titanium, square head	4
V00050091A	Screw Ø 1.5 x 18 locked, titanium, square head	4
V00050092A	Screw Ø 1.5 x 20 locked, titanium, square head	4
V00040129A	Screw Ø 2.0 x 6 cortical, titanium, square head	4
V00040103A	Screw Ø 2.0 x 8 cortical, titanium, square head	4
V00040104A	Screw Ø 2.0 x 10 cortical, titanium, square head	4
V00040105A	Screw Ø 2.0 x 12 cortical, titanium, square head	4
V00040106A	Screw Ø 2.0 x 14 cortical, titanium, square head	4
V00040107A	Screw Ø 2.0 x 16 cortical, titanium, square head	4
V00040133A	Screw Ø 2.0 x 18 cortical, titanium, square head	4
V00040134A	Screw Ø 2.0 x 20 cortical, titanium, square head	4
V00040130A	Screw Ø 1.5 x 6 cortical, titanium, square head	4
V00040108A	Screw Ø 1.5 x 8 cortical, titanium, square head	4
V00040109A	Screw Ø 1.5 x 10 cortical, titanium, square head	4
V00040110A	Screw Ø 1.5 x 12 cortical, titanium, square head	4
V00040111A	Screw Ø 1.5 x 14 cortical, titanium, square head	4
V00040112A	Screw Ø 1.5 x 16 cortical, titanium, square head	4
V00040132A	Screw Ø 1.5 x 16 cortical, titanium, square head	4
V00040133A	Screw Ø 1.5 x 16 cortical, titanium, square head	4

Osti-Lok locked plate system

OSTI-LOK MINI

Osti-Lok mini set

Osti-Lok mini set



Suited for the treatment of patients up to 10 kg

€ 5,270.54

Kit composition

Code	Product description	Quantity
U00010045A	Quick coupling with 1.0 mm K-wire	1
U00010044A	Quick coupling with 1.5 mm K-wire	1
U00010046A	Mini Screwdriver for square head screws	1
M00070003A	Mini depth gauge	1
S00050005A	Drill guide, mini	2
P00040013A	Angled grip-nut, mini	1
P00040014A	Straight grip nut, mini	2
U00010047A	Grip nut handle, mini	2
U00010048A	In-plane bender, mini	2
U00010049A	Hex-key	1
U00010052A	Deiver insert, mini	1
C00040031A	Mini Osti-lok tray	1

Osti-Lok locked plate system

OSTI-LOK MINI

Locking plates Osti-Lok mini

**Straight plate,
thickness 2.0,
19 holes**



I00020071A

Titanium

€ 201.30

**Straight plate,
thickness 1.5,
19 holes**



I00020074A

Titanium

€ 192.50

**Y plate,
thickness 2.0,
8 holes**



I00020085A

Titanium

€ 183.70

**Y plate,
thickness 1.5,
8 holes**



I00020087A

Titanium

€ 180.40

**Cross plate,
thickness 2.0,
13 holes**



I00020069A

Titanium

€ 183.70

**Cross plate,
thickness 1.5,
13 holes**



I00020072A

Titanium

€ 180.40

**High strength plate,
thickness 2.0,
10 holes**



I00020084A

Titanium

€ 154.00

**High strength plate,
thickness 1.5,
10 holes**

I00020086A

Tiatanio

€ 150.70

Osti-Lok locked plate system

OSTI-LOK MINI

Locking screws Osti-Lok mini

**Ø 2.0 x 6 locking screw,
titanium,
square head**

V00050073A

Titanium

€ 23.76

**Ø 2.0 x 8 locking screw,
titanium,
square head**



V00050038A

Titanium

€ 23.76

**Ø 2.0 x 10 locking screw,
titanium,
square head**

V00050039A

Titanium

€ 23.76

**Ø 2.0 x 12 locking screw,
titanium,
square head**

V00050040A

Titanium

€ 23.76

**Ø 2.0 x 14 locking screw,
titanium,
square head**

V00050041A

Titanium

€ 23.76

**Ø 2.0 x 16 locking screw,
titanium,
square head**

V00050042A

Titanium

€ 23.76

**Ø 2.0 x 18 locking screw,
titanium,
square head**

V00050093A

Titanium

€ 23.76

**Ø 2.0 x 20 locking screw,
titanium,
square head**

V00050094A

Titanium

€ 23.76

**Ø 1.5 x 6 locking screw,
titanium,
square head**

V00050074A

Titanium

€ 23.76

**Ø 1.5 x 8 locking screw,
titanium,
square head**

V00050043A

Titanium

€ 23.76

**Ø 1.5 x 10 locking screw,
titanium,
square head**

V00050044A

Titanium

€ 23.76



**Ø 1.5 x 12 locking screw,
titanium,
square head**

V00050045A

Titanium

€ 23.76

**Ø 1.5 x 14 locking screw,
titanium,
square head**

V00050046A

Titanium

€ 23.76

**Ø 1.5 x 16 locking screw,
titanium,
square head**



V00050047A

Titanium

€ 23.76

**Ø 1.5 x 18 locking screw,
titanium,
square head**

V00050091A

Titanium

€ 23.76

**Ø 1.5 x 20 locking screw,
titanium,
square head**

V00050092A

Titanium

€ 23.76

Osti-Lok locked plate system

OSTI-LOK MINI

Cortical screws Osti-Lok mini

**Ø 2.0 x 6 cortical screw,
titanium,
square head**

V00040129A

Titanium

€ 20.14

**Ø 2.0 x 8 cortical screw,
titanium,
square head**

V00040103A

Titanium

€ 20.14

**Ø 2.0 x 10 cortical screw,
titanium,
square head**



V00040104A

Titanium

€ 20.14

**Ø 2.0 x 12 cortical screw,
titanium,
square head**

V00040105A

Titanium

€ 20.14

**Ø 2.0 x 14 cortical screw,
titanium,
square head**

V00040106A

Titanium

€ 20.14

**Ø 2.0 x 16 cortical screw,
titanium,
square head**

V00040107A

Titanium

€ 20.14

**Ø 2.0 x 18 cortical screw,
titanium,
square head**



V00040133A

Titanium

€ 20.14

**Ø 2.0 x 20 cortical screw,
titanium,
square head**

V00040134A

Titanium

€ 20.14

**Ø 1.5 x 6 cortical screw,
titanium,
square head**

V00040130A

Titanium

€ 20.14

**Ø 1.5 x 8 cortical screw,
titanium,
square head**

V00040108A

Titanium

€ 20.14

**Ø 1.5 x 10 cortical screw,
titanium,
square head**



V00040109A

Titanium

€ 20.14

**Ø 1.5 x 12 cortical screw,
titanium,
square head**

V00040110A

Titanium

€ 20.14

**Ø 1.5 x 14 cortical screw,
titanium,
square head**

V00040111A

Titanium

€ 20.14

Ø 1.5 x 16 cortical screw,
titanium,
square head



V00040112A

Titanium

€ 20.14

Ø 1.5 x 18 cortical screw,
titanium,
square head




V00040131A

Titanio

€ 20.14

Ø 1.5 x 20 cortical screw,
titanium,
square head



V00040132A

Titanio

€ 20.14

Osti-Lok locked plate system

OSTI-LOK MINI

Osti-Lok Mini replacement elements

Square driver, mini



U00010046A

Titanium

€ 99.00

Deiver insert, mini



U00010052A

Titanium

€ 44.00

Depth gauge, mini



M00070003A

Titanium

€ 110.00

Angled grip-nut, mini



P00040013A

Titanium

€ 44.00

Grip nut, mini



P00040014A

Titanium

€ 44.00

Grip nut handle, mini



U00010047A

Titanium

€ 44.00

AO coupling for 1,5 mm K-wire



U00010044A

Titanium

€ 44.00

AO coupling for 1,0 mm K-wire



U00010045A

Titanium

€ 44.00

In-plane bender, mini



U00010048A

Titanium

€ 44.00

Hex-key, 3/32



U00010049A

Titanium

€ 1.10

**Drill guide, mini
with lock**



S00050006A

Titanium

€ 44.00

Mini Osti-Lok tray



C00040031A

Titanium

€ 490.80

AD  **MAIORA**®
WE MOVE ON



External Fixation

UNIVERSAL SYSTEM

Mini linear radiolucent fixator kit

Polilock Fixator

**Mini linear radiolucent
fixator kit
Polilock Fixator**



K00010001B

*Mini radiolucent set for the
treatment of patients
weighing up to 5 kg*

€ 316.50

Kit composition

Code	Product description	Quantity
B00010001A	Carbon rod Ø 5.0 mm Length 100 mm	2
B00010002A	Carbon rod Ø 5.0 mm Length 150 mm	2
B00010010A	Carbon rod Ø 5.0 mm Length 50 mm	2
B00010011A	Carbon rod Ø 5.0 mm Length 75 mm	2
B00010012A	Carbon rod Ø 5.0 mm Length 125 mm	2
C00040008A	Autoclavable complete container for linear radiolucent mini kit	1
D00010019A	Nut M3 ss self-locking	20
F00020173A	Threaded pin Ø 2.0 mm Length 100 mm, fil. 15 mm	2
F00020174A	Threaded pin Ø 2.4 mm Length 100 mm, fil. 15 mm	2
F00020182A	K wire Ø 1,2 mm Length 120 mm, negative 10 mm thread	2
F00020183A	K wire Ø 1,6 mm Length 120 mm, negative 10 mm thread	4
G00010001A	Hook Ø 3 mm Length 16.5 mm x 5.0 mm	10
G00010002A	Hook Ø 3 mm Length 24.5 mm x 5.0 mm	6
G00010004A	Hook Ø 3 mm Length 16.5 mm x 7.5 mm	6
M00010020A	Round clamp Ø 20 mm 3H 3 mm aligned B5 rxl black	8
M00010021A	Round clamp Ø 20 mm 3H 3 mm staggered B5 rxl black	4
P00040004A	Cap for pins Ø 1.6 mm	8
P00040005A	Cap for pins Ø 1.8 mm	4
R00020006A	Washer ss toothed hole 3 mm	20
V00020013A	Screw hexagonal head ss M3x20 mm	10

UNIVERSAL SYSTEM

Rail fixator Kit

Rail fixator Kit



K00010014A

*Suitable for patients from
10 to 40 kg in weight*

€ 451.04


Kit composition

Code	Product description	Quantity
B00040011A	Rail 2 holes radiolucent	1
B00040012A	Rail 3 holes radiolucent	1
B00040013A	Rail 4 holes radiolucent	1
C00040025A	Autoclavable complete container for rail fixator kit	1
D00010018A	M6 nut for clamp and post	8
M00010003A	Wire/pin clamp aluminum alloy	4
M00010026A	Double-purchase post, M6	4
R00020003A	Round washer 15x5.5 mm	8
R00030006A	Rectangular washer 8.5x15 mm	4
R00030009A	Rectangular self-centering steel washer 8.2x15 mm	4
S00040010A	Distractor 16 cm	1

UNIVERSAL SYSTEM

Dynamizable kit

Dynamizable kit



K00010002A

Suitable for patients weighing more than 10 kg

€ 553.15

Kit composition

Code	Product description	Quantity
C00040010A	Autoclavable complete container for dynamizable fixator kit	1
E00010001A	Clamp extension for dynamizable	2
F00020175A	Threaded pin Ø 2.7 mm Length 120 mm, fil. 20 mm	4
F00020176A	Threaded pin Ø 3.5 mm Length 130 mm, fil. 25 mm	4
F00030001B	Dynamizable linear fixator 10 cm	1
F00030002B	Dynamizable linear fixator 15 cm	1
M00010018A	Lightweight clamp for dynamizable fixator	6
P00040002A	Cap for pins Ø 3,0 mm	4
P00040005A	Cap for pins Ø 1,8 mm	4
U00010002A	Hexagon wrench 3.0 mm with spherical tip	1

UNIVERSAL SYSTEM

Mini hybrid radiolucent fixator kit

Midi hybrid radiolucent fixator kit



K00010003A

*Recommended use:
small dogs and cats*

€1,050.30

Kit composition

Code	Product description	Quantity
A00020005B	Carbon ring 270° Ø 54 mm x 3 mm	2
A00020009B	Carbon ring 360° Ø 54 mm x 3 mm	3
B00020015A	Aluminum M5 threaded rod Length 50 mm	4
B00020019A	Aluminum M5 threaded rod Length 130 mm	4
B00020035A	Aluminum M5 threaded rod Length 70 mm	4
B00020036A	Aluminum M5 threaded rod Length 90 mm	4
B00020037A	Aluminum M5 threaded rod Length 110 mm	4
B00040002A	Rail 1 oval hole, inclinable, radiolucent, M6 thread	1
B00040004A	Rail 1 oval hole, orthogonal, radiolucent, M6 thread	1
C00020001A	Pin-holder cylinder 1 hole	1
C00020002A	Pin-holder cylinder 2 holes	1
C00020003A	Pin-holder cylinder 3 holes	1
C00040018A	Complete autoclavable container for mini hybrid kit	1
D00010003A	M5 small nut with flange	20
D00010008A	Nut M6 aluminum gold self-locking	1
D00010018A	M6 nut for clamp and post	3
D00010019A	Nut M3 ss selflocking	20
D00020003A	Hemispherical nut M5 aluminum	12
D00020004A	M6 hemispherical nut aluminum alloy	2
F00020157A	K wire half tip, Ø 1.0 mm Length 220 mm	5
F00020158A	K wire half tip, Ø 1.2 mm Length 250 mm	5
F00020173A	Threaded pin Ø 2.0 mm Length 100 mm, fil. 15 mm	4
F00020174A	Threaded pin Ø 2.4 mm Length 100 mm, fil. 15 mm	4
F00020183A	K wire Ø 1,6 mm Length 120 mm, negative 10 mm thread	2
G00010001A	Hook Ø 3 mm Length 16.5 mm x 5.0 mm	4
G00010002A	Hook Ø 3 mm Length 24.5 mm x 5.0 mm	4
G00010004A	Hook Ø 3 mm Length 16.5 mm x 7.5 mm	4
M00010003A	Wire/pin clamp aluminum alloy	3
M00010020A	Round clamp Ø 20 mm 3H 3 mm aligned B5 rxl black	4
M00010021A	Round clamp Ø 20 mm 3H 3 mm staggered B5 rxl black	4
P00040004A	Cap for pins Ø 1.6 mm	2
P00040005A	Cap for pins Ø 1,8 mm	8
R00020002A	Round washer 14x7.5 mm	3
R00020003A	Round washer 15x5.5 mm	20
R00030001A	Hemispheric washer inox 10x16 mm	2
R00030002A	Hemispheric washer inox 8x14 mm	12
R00030004A	Hemispherical washer rxlucent 10x15	12
R00030008A	Rectangular self-centering steel washer 5.5x15 mm	24
R00030009A	Rectangular self-centering steel washer 8.2x15 mm	6
V00020001A	Threaded locker for pin Length 8 mm	4
V00020002A	Threaded locker for pin Length 12 mm	4
V00020003A	Threaded locker for pin Length 16 mm	4
V00030002A	Hexagonal head perforated screw M5 Length 13,5 mm, for K wires	12

UNIVERSAL SYSTEM

Mini hybrid aluminum fixator kit

Mini hybrid aluminum external fixator kit



K00010013A

Suitable for the treatment of patients up to 15 kg

€ 963.18

Kit composition

Code	Product description	Quantity
A00010004A	Aluminum alloy ring 270° Ø 55 mm x 6 mm	2
A00010007A	Aluminum alloy ring 360° Ø 55 mm x 6 mm	3
B00020015A	Aluminum M5 threaded rod Length 50 mm	4
B00020019A	Aluminum M5 threaded rod Length 130 mm	4
B00020035A	Aluminum M5 threaded rod Length 70 mm	4
B00020036A	Aluminum M5 threaded rod Length 90 mm	4
B00020037A	Aluminum M5 threaded rod Length 110 mm	4
B00030001A	Rail with 1 oval hole, inclined positioning, aluminum alloy	1
B00030003A	Rail with 1 oval hole, orthogonal, aluminum alloy	1
C00020001A	Pin-holder cylinder 1 hole	1
C00020002A	Pin-holder cylinder 2 holes	1
C00020003A	Pin-holder cylinder 3 holes	1
C00040018A	Complete autoclavable container for mini hybrid kit	1
D00010003A	M5 small nut with flange	20
D00010008A	Nut M6 aluminum gold self-locking	1
D00010018A	M6 nut for clamp and post	3
D00010019A	Nut M3 ss selflocking	20
D00020003A	Hemispherical nut M5 aluminum	12
D00020004A	M6 hemispherical nut aluminum alloy	2
F00020157A	K wire half tip, Ø 1.0 mm Length 220 mm	5
F00020158A	K wire half tip, Ø 1.2 mm Length 250 mm	5
F00020173A	Threaded pin Ø 2.0 mm Length 100 mm, fil. 15 mm	4
F00020174A	Threaded pin Ø 2.4 mm Length 100 mm, fil. 15 mm	4
F00020183A	K wire Ø 1,6 mm Length 120 mm, negative 10 mm thread	2
G00010001A	Hook Ø 3 mm Length 16.5 mm x 5.0 mm	4
G00010002A	Hook Ø 3 mm Length 24.5 mm x 5.0 mm	4
G00010004A	Hook Ø 3 mm Length 16.5 mm x 7.5 mm	4
M00010003A	Wire/pin clamp aluminum alloy	3
M00010020A	Round clamp Ø 20 mm 3H 3 mm aligned B5 rxl black	4
M00010021A	Round clamp Ø 20 mm 3H 3 mm staggered B5 rxl black	4
P00040004A	Cap for pins Ø 1.6 mm	2
P00040005A	Cap for pins Ø 1,8 mm	8
R00020002A	Round washer 14x7.5 mm	3
R00020003A	Round washer 15x5.5 mm	20
R00030001A	Hemispheric washer inox 10x16 mm	2
R00030002A	Hemispheric washer inox 8x14 mm	12
R00030003A	Hemispherical washer rxlucent 10x15	12
R00030008A	Rectangular self-centering steel washer 5.5x15 mm	24
R00030009A	Rectangular self-centering steel washer 8.2x15 mm	6
V00020001A	Threaded locker for pin Length 8 mm	4
V00020002A	Threaded locker for pin Length 12 mm	4
V00020003A	Threaded locker for pin Length 16 mm	4
V00030002A	Hexagonal head perforated screw M5 Length 13,5 mm, for K wires	12

UNIVERSAL SYSTEM

Midi hybrid radiolucent fixator kit

Midi hybrid radiolucent fixator kit



K00010005A

Recommended use:
dogs from 10 to 35 kg

€1,367.50

Kit composition

Code	Product description	Quantity
A00020007B	Carbon ring 270° Ø 84 mm x 6 mm	2
A00020011B	Carbon ring 360° Ø 84 mm x 6 mm	3
B00020032A	Threaded rod M5 Length 170 mm, steel	4
B00020014A	Threaded rod M5 Length 190 mm, steel	4
B00020026A	Threaded rod M5 Length 50 mm, steel	4
B00020027A	Threaded rod M5 Length 70 mm, steel	4
B00020028A	Threaded rod M5 Length 90 mm, steel	4
B00020029A	Threaded rod M5 Length 110 mm, steel	4
B00020030A	Threaded rod M5 Length 130 mm, steel	4
B00020031A	Threaded rod M5 Length 150 mm, steel	4
B00040002A	Rail 1 oval hole, inclinable, radiolucent, M6 thread	1
B00040004A	Rail 1 oval hole, orthogonal, radiolucent, M6 thread	1
B00040006A	Rail 2 oval holes, inclinable, radiolucent, M6 thread	1
B00040008A	Rail 2 oval holes, orthogonal, radiolucent, M6 thread	1
C00020001A	Pin-holder cylinder 1 hole	1
C00020002A	Pin-holder cylinder 2 holes	1
C00020003A	Pin-holder cylinder 3 holes	1
C00040016A	Complete autoclavable container for linear radiolucent midi kit	1
D00010003A	M5 small nut with flange	20
D00010008A	Nut M6 aluminum gold self-locking	1
D00010018A	M6 nut for clamp and post	3
D00010019A	Nut M3 ss selflocking	20
D00020003A	Hemispherical nut M5 aluminum	12
D00020004A	M6 hemispherical nut aluminum alloy	2
F00020158A	K wire half tip, Ø 1.2 mm Length 250 mm	6
F00020159A	K wire half tip, Ø 1.5 mm Length 250 mm	6
F00020174A	Threaded pin Ø 2.4 mm Length 100 mm, fil. 15 mm	4
F00020175A	Threaded pin Ø 2.7 mm Length 120 mm, fil. 20 mm	4
F00020176A	Threaded pin Ø 3.5 mm Length 130 mm, fil. 25 mm	4
G00010001A	Hook Ø 3 mm Length 16.5 mm x 5.0 mm	4
G00010004A	Hook Ø 3 mm Length 16.5 mm x 7.5 mm	4
M00010003A	Wire/pin clamp aluminum alloy	3
M00010020A	Round clamp Ø 20 mm 3H 3 mm aligned B5 rxl black	4
M00010021A	Round clamp Ø 20 mm 3H 3 mm staggered B5 rxl black	4
P00040002A	Cap for pins Ø 3,0 mm	4
P00040005A	Cap for pins Ø 1,8 mm	8
R00020002A	Round washer 14x7.5 mm	3
R00020003A	Round washer 15x5.5 mm	20
R00030001A	Hemispheric washer inox 10x16 mm	4

UNIVERSAL SYSTEM

Midi hybrid radiolucent fixator kit

Midi hybrid radiolucent fixator kit



K00010005A

*Recommended use:
dogs from 10 to 35 kg*

€ 1,367.50

Kit composition

Code	Product description	Quantity
R00030002A	Hemispheric washer inox 8x14 mm	12
R00030008A	Rectangular self-centering steel washer 5.5x15 mm	24
R00030009A	Rectangular self-centering steel washer 8.2x15 mm	6
V00020001A	Threaded locker for pin Length 8 mm	4
V00020002A	Threaded locker for pin Length 12 mm	4
V00020003A	Threaded locker for pin Length 16 mm	4
V00030001A	Hexagonal head perforated screw M5 Length 13,5 mm, for K wires	12

UNIVERSAL SYSTEM

Midi hybrid aluminum fixator kit

Midi hybrid aluminum fixator kit



K00010006A

*Recommended use:
Dogs from 10 to 40 kg*

€1,087.52

Kit composition

Code	Product description	Quantity
A00010005A	Aluminum alloy ring 270° Ø 85 mm x 6 mm	2
A00010008A	Aluminum alloy ring 360° Ø 85 mm x 6 mm	3
B00020032A	Threaded rod M5 Length 170 mm, steel	4
B00020014A	Threaded rod M5 Length 190 mm, steel	4
B00020026A	Threaded rod M5 Length 50 mm, steel	4
B00020027A	Threaded rod M5 Length 70 mm, steel	4
B00020028A	Threaded rod M5 Length 90 mm, steel	4
B00020029A	Threaded rod M5 Length 110 mm, steel	4
B00020030A	Threaded rod M5 Length 130 mm, steel	4
B00020031A	Threaded rod M5 Length 150 mm, steel	4
B00030001A	Rail with 1 oval hole, inclined positioning, aluminum alloy	1
B00030003A	Rail with 1 oval hole, orthogonal, aluminum alloy	1
B00030004A	Rail with 2 oval holes, inclined positioning, aluminum alloy	1
B00030005A	Rail with 2 oval holes, orthogonal, aluminum alloy	1
C00020001A	Pin-holder cylinder 1 hole	1
C00020002A	Pin-holder cylinder 2 holes	1
C00020003A	Pin-holder cylinder 3 holes	1
C00040013A	Complete autoclavable container for linear midi aluminum kit	1
D00010003A	M5 small nut with flange	20
D00010008A	Nut M6 aluminum gold self-locking	1
D00010018A	M6 nut for clamp and post	3
D00010019A	Nut M3 ss selflocking	20
D00020003A	Hemispherical nut M5 aluminum	12
D00020006A	M7 hemispherical nut aluminum alloy	2
F00020158A	K wire half tip, Ø 1.2 mm Length 250 mm	6
F00020159A	K wire half tip, Ø 1.5 mm Length 250 mm	6
F00020174A	Threaded pin Ø 2.4 mm Length 100 mm, fil. 15 mm	4
F00020175A	Threaded pin Ø 2.7 mm Length 120 mm, fil. 20 mm	4
F00020176A	Threaded pin Ø 3.5 mm Length 130 mm, fil. 25 mm	4
G00010001A	Hook Ø 3 mm Length 16.5 mm x 5.0 mm	4
G00010004A	Hook Ø 3 mm Length 16.5 mm x 7.5 mm	4
M00010003A	Wire/pin clamp aluminum alloy	3
M00010020A	Round clamp Ø 20 mm 3H 3 mm aligned B5 rxl black	4
M00010021A	Round clamp Ø 20 mm 3H 3 mm staggered B5 rxl black	4
P00040002A	Cap for pins Ø 3,0 mm	4
P00040005A	Cap for pins Ø 1,8 mm	8
R00020002A	Round washer 14x7.5 mm	3
R00020003A	Round washer 15x5.5 mm	20
R00030001A	Hemispheric washer inox 10x16 mm	4

External Fixation

UNIVERSAL SYSTEM

Midi hybrid aluminum fixator kit

Midi hybrid aluminum fixator kit



K00010006A

*Recommended use:
Dogs from 10 to 40 kg*

€ 1,087.52

Kit composition

Code	Product description	Quantity
R00030002A	Hemispheric washer inox 8x14 mm	12
R00030008A	Rectangular self-centering steel washer 5.5x15 mm	24
R00030009A	Rectangular self-centering steel washer 8.2x15 mm	6
V00020001A	Threaded locker for pin Length 8 mm	4
V00020002A	Threaded locker for pin Length 12 mm	4
V00020003A	Threaded locker for pin Length 16 mm	4
V00030001A	Hexagonal head perforated screw M5 Length 13,5 mm, for K wires	12

UNIVERSAL SYSTEM

Maxi hybrid radiolucent fixator kit

Maxi hybrid radiolucent fixator kit



K00010015A

Suitable for the treatment of patients from 15 to 50 kg

€1,488.54

Kit composition

Code	Product description	Quantity
A00020016A	Carbon ring 360° Ø 114 mm x 9 mm	3
B00020032A	Threaded rod M5 Length 170 mm, steel	4
B00020014A	Threaded rod M5 Length 190 mm, steel	4
B00020026A	Threaded rod M5 Length 50 mm, steel	4
B00020027A	Threaded rod M5 Length 70 mm, steel	4
B00020028A	Threaded rod M5 Length 90 mm, steel	4
B00020029A	Threaded rod M5 Length 110 mm, steel	4
B00020030A	Threaded rod M5 Length 130 mm, steel	4
B00020031A	Threaded rod M5 Length 150 mm, steel	4
B00040002A	Rail 1 oval hole, inclinable, radiolucent, M6 thread	1
B00040004A	Rail 1 oval hole, orthogonal, radiolucent, M6 thread	1
B00040006A	Rail 2 oval holes, inclinable, radiolucent, M6 thread	1
B00040008A	Rail 2 oval holes, orthogonal, radiolucent, M6 thread	1
B00040009A	Rail 3 oval hole, inclinable, radiolucent, M6 thread	1
B00040010A	Rail 3 oval hole, orthogonal, radiolucent, M6 thread	1
C00020001A	Pin-holder cylinder 1 hole	1
C00020002A	Pin-holder cylinder 2 holes	1
C00020003A	Pin-holder cylinder 3 holes	1
C00040015A	Complete autoclavable container for kit hybrid maxi	1
D00010003A	M5 small nut with flange	20
D00010018A	M6 nut for clamp and post	3
D00010019A	Nut M3 ss selflocking	20
D00020003A	Hemispherical nut M5 aluminum	12
D00020006A	M7 hemispherical nut aluminum alloy	2
F00020158A	K wire half tip, Ø 1.2 mm Length 250 mm	6
F00020159A	K wire half tip, Ø 1.5 mm Length 250 mm	6
F00020175A	Threaded pin Ø 2.7 mm Length 120 mm, fil. 20 mm	4
F00020176A	Threaded pin Ø 3.5 mm Length 130 mm, fil. 25 mm	4
F00020177A	Threaded pin Ø 4,5 mm Length 130 mm, fil. 25 mm	4
G00010001A	Hook Ø 3 mm Length 16.5 mm x 5.0 mm	4
G00010004A	Hook Ø 3 mm Length 16.5 mm x 7.5 mm	4
M00010003A	Wire/pin clamp aluminum alloy	3
M00010020A	Round clamp Ø 20 mm 3H 3 mm aligned B5 rxl black	4
M00010021A	Round clamp Ø 20 mm 3H 3 mm staggered B5 rxl black	4
P00040002A	Cap for pins Ø 3,0 mm	4
P00040003A	Cap for pins Ø 4,0 mm	4
P00040005A	Cap for pins Ø 1,8 mm	8
R00020002A	Round washer 14x7.5 mm	3
R00020003A	Round washer 15x5.5 mm	20

UNIVERSAL SYSTEM

Maxi hybrid radiolucent fixator kit

Maxi hybrid radiolucent fixator kit



K00010015A

Suitable for the treatment of patients from 15 to 50 kg

€ 1,488.54

Kit composition

Code	Product description	Quantity
R00030001A	Hemispheric washer inox 10x16 mm	4
R00030002A	Hemispheric washer inox 8x14 mm	12
R00030008A	Rectangular self-centering steel washer 5.5x15 mm	24
R00030009A	Rectangular self-centering steel washer 8.2x15 mm	6
V00020001A	Threaded locker for pin Length 8 mm	4
V00020002A	Threaded locker for pin Length 12 mm	4
V00020003A	Threaded locker for pin Length 16 mm	4
V00030001A	Hexagonal head perforated screw M5 Length 13,5 mm, for K wires	12

UNIVERSAL SYSTEM

Maxi hybrid aluminum fixator kit

Maxi hybrid aluminum fixator kit



K00010007A

Indicated for the treatment of patients from 15 to 50 kg

€1.180,04

Kit composition

Code	Product description	Quantity
A00010006A	Aluminum alloy ring 270° Ø 115 mm x 6 mm	2
A00010009A	Aluminum alloy ring 360° Ø 115 mm x 6 mm	3
B00020032A	Threaded rod M5 Length 170 mm, steel	4
B00020014A	Threaded rod M5 Length 190 mm, steel	4
B00020026A	Threaded rod M5 Length 50 mm, steel	4
B00020027A	Threaded rod M5 Length 70 mm, steel	4
B00020028A	Threaded rod M5 Length 90 mm, steel	4
B00020029A	Threaded rod M5 Length 110 mm, steel	4
B00020030A	Threaded rod M5 Length 130 mm, steel	4
B00020031A	Threaded rod M5 Length 150 mm, steel	4
B00030001A	Rail with 1 oval hole, inclined positioning, aluminum alloy	1
B00030003A	Rail with 1 oval hole, orthogonal, aluminum alloy	1
B00030004A	Rail with 2 oval holes, inclined positioning, aluminum alloy	1
B00030005A	Rail with 2 oval holes, orthogonal, aluminum alloy	1
B00030006A	Rail with 3 oval hole, inclined positioning, aluminum alloy	1
B00030007A	Rail with 3 oval hole, orthogonal, aluminum alloy	1
C00020001A	Pin-holder cylinder 1 hole	1
C00020002A	Pin-holder cylinder 2 holes	1
C00020003A	Pin-holder cylinder 3 holes	1
C00040015A	Complete autoclavable container for kit hybrid maxi	1
D00010003A	M5 small nut with flange	20
D00010018A	M6 nut for clamp and post	3
D00010019A	Nut M3 ss selflocking	20
D00020003A	Hemispherical nut M5 aluminum	12
D00020006A	M7 hemispherical nut aluminum alloy	2
F00020158A	K wire half tip, Ø 1.2 mm Length 250 mm	6
F00020159A	K wire half tip, Ø 1.5 mm Length 250 mm	6
F00020175A	Threaded pin Ø 2.7 mm Length 120 mm, fil. 20 mm	4
F00020176A	Threaded pin Ø 3.5 mm Length 130 mm, fil. 25 mm	4
F00020177A	Threaded pin Ø 4,5 mm Length 130 mm, fil. 25 mm	4
G00010001A	Hook Ø 3 mm Length 16.5 mm x 5.0 mm	4
G00010004A	Hook Ø 3 mm Length 16.5 mm x 7.5 mm	4
M00010003A	Wire/pin clamp aluminum alloy	3
M00010020A	Round clamp Ø 20 mm 3H 3 mm aligned B5 rxl black	4
M00010021A	Round clamp Ø 20 mm 3H 3 mm staggered B5 rxl black	4
P00040002A	Cap for pins Ø 3,0 mm	4
P00040003A	Cap for pins Ø 4,0 mm	4
P00040005A	Cap for pins Ø 1,8 mm	8
R00020002A	Round washer 14x7.5 mm	3
R00020003A	Round washer 15x5.5 mm	20

UNIVERSAL SYSTEM

Maxi hybrid aluminum fixator kit

Maxi hybrid aluminum fixator kit



K00010007A

Indicated for the treatment of patients from 15 to 50 kg

€ 1.180,04

Kit composition

Code	Product description	Quantity
R00030001A	Hemispheric washer inox 10x16 mm	4
R00030002A	Hemispheric washer inox 8x14 mm	12
R00030008A	Rectangular self-centering steel washer 5.5x15 mm	24
R00030009A	Rectangular self-centering steel washer 8.2x15 mm	6
V00020001A	Threaded locker for pin Length 8 mm	4
V00020002A	Threaded locker for pin Length 12 mm	4
V00020003A	Threaded locker for pin Length 16 mm	4
V00030001A	Hexagonal head perforated screw M5 Length 13,5 mm, for K wires	12

UNIVERSAL SYSTEM

Rings

Aluminum alloy ring
360° Ø 115 mm x 6 mm



A0001009A

Aluminum

€ 36.30

Aluminum alloy ring
360° Ø 85 mm x 6 mm



A0001008A

Aluminum

€ 33.00

Aluminum alloy ring
360° Ø 55 mm x 6 mm



A0001007A

Aluminum

€ 27.50

Aluminum alloy ring
270° Ø 115 mm x 6 mm



A0001006A

Aluminum

€ 34.10

Aluminum alloy ring
270° Ø 85 mm x 6 mm



A0001005A

Aluminum

€ 30.80

Aluminum alloy ring
270° Ø 55 mm x 6 mm



A0001004A

Aluminum

€ 27.50

Aluminum alloy ring
180° Ø 115 mm x 6 mm



A0001003A

Aluminum

€ 25.30

Aluminum alloy ring
180° Ø 85 mm x 6 mm



A0001002A

Aluminum

€ 22.00

Aluminum alloy ring
180° Ø 55 mm x 6 mm



A0001001A

Aluminum

€ 18.70

Carbon ring
360° Ø 114 mm x 9 mm



A0002016A

Carbon composite

€ 165.00

Carbon ring
360° Ø 84 mm x 6 mm

A0002011B

Carbon composite

€ 93.50

Carbon ring
360° Ø 54 mm x 3 mm

A0002009B

Carbon composite

€ 46.20

Carbon ring
270° Ø 84 mm x 6 mm



A0002007B

Carbon composite

€ 88.00

Carbon ring
270° Ø 54 mm x 3 mm



A0002005B

Carbon composite

€ 44.00

External Fixation

UNIVERSAL SYSTEM

Threaded bars

**Aluminum M5
threaded bar**
Length 50 mm

B00020034A

Aluminum

€ 4.62

**Aluminum M5
threaded bar**
Length 70 mm



B00020035A

Aluminum

€ 4.62

**Aluminum M5
threaded bar**
Length 90 mm



B00020036A

Aluminum

€ 4.62

**Aluminum M5
threaded bar**
Length 110 mm

B00020037A

Aluminum

€ 4.95

**Aluminum M5
threaded bar**
Length 130 mm

B00020019A

Aluminum

€ 4.95

**Stainless steel
threaded bar M5
Length 50 mm**

B00020026A

Stainless steel

€ 3.30

**Stainless steel
threaded bar M5
Length 70 mm**

B00020027A

Stainless steel

€ 3.52

**Stainless steel
threaded bar M5
Length 90 mm**

B00020028A

Stainless steel

€ 3.96

**Stainless steel
threaded bar M5
Length 110 mm**

B00020029A

Stainless steel

€ 4.40

**Stainless steel
threaded bar M5
Length 130 mm**

B00020030A

Stainless steel

€ 4.62

**Stainless steel
threaded bar M5
Length 150 mm**

B00020031A

Stainless steel

€ 4.84

**Stainless steel
threaded bar M5
Length 170 mm**

B00020032A

Stainless steel

€ 5.28



**Stainless steel
threaded bar M5
Length 190 mm**



B00020033A

Stainless steel

€ 5.50

External Fixation

UNIVERSAL SYSTEM Rails

**Rail 1 oval hole,
inclined, radiolucent,
M6 thread**



B00040002A

Technopolymer

€ 39.60

**Rail 1 oval hole,
orthogonal, radiolucent,
M6 thread**



B00040004A

Technopolymer

€ 37.40

**Rail with 1 oval hole,
inclined**



B00030001A

Aluminum

€ 35.64

**Rail with 1 oval hole,
orthogonal**



B00030003A

Aluminum

€ 33.66

**Rail 2 oval holes,
inclined, radiolucent,
M6 thread**



B00040006A

Technopolymer

€ 44.00

**Rail 2 oval holes,
orthogonal,
radiolucent, M6 thread**



B00040008A

Technopolymer

€ 44.00

**Rail 2 holes radiolucent
for rail fixator**



B00040011A

Technopolymer

€ 38.50

**Rail with 2 oval holes,
inclinable aluminum
alloy**



B00030004A

Aluminum

€ 40.26

**Rail with 2 oval holes,
orthogonal, aluminum
alloy**



B00030005A

Aluminum

€ 40.26

**Rail 3 oval holes,
inclinable, radiolucent,
M6 thread**



B00040009A

Technopolymer

€ 48.40

**Rail 3 oval holes,
orthogonal,
radiolucent, M6 thread**



B00040010A

Technopolymer

€ 48.40

**Rail 3 holes radiolucent
for rail fixator**



B00040012A

Technopolymer

€ 44.00

**Rail 1 oval hole,
inclined, radiolucent,
M6 thread**



B00030007A

Aluminum

€ 43.56

**Rail with 3 oval holes,
orthogonal, aluminum
alloy**



B00030006A

Aluminum

€ 43.56

**Rail 4 holes
radiolucent for rail
fixator**



B00040013A

Technopolymer

€ 49.50

**Inclinable cylindrical rod
Length 80 mm
For M00010018A clamps**



A00030001A

Stainless steel

€ 26.40

**Inclinable cylindrical rod
Length 110 mm
For M00010018A clamps**



A00030002A

Stainless steel

€ 30.80

UNIVERSAL SYSTEM

Nuts

**Nut M5 silver
self-locking aluminum**



D00010007A

Aluminum

€ 1.49

Nut M3 self-locking



D00010019A

Stainless steel

€ 0.33

M3 nut



D00010010A

Stainless steel

€ 0.17

**M6 nut for clamp and
post**



D00010018A

Stainless steel

€ 0.55

**Nut M6 aluminum gold
self-locking**



D00010008A

Aluminum

€ 1.49

M4 hemispherical nut



D00020001A

Aluminum

€ 8.80

M5 hemispherical nut

D00020003A

Aluminum

€ 8.80

M5 hemispherical nut



D00020002A

Stainless steel

€ 13.20

M5 hemispherical nut

D00020010A

Titanium

€ 19.80

M6 hemispherical nut



D00020004A

Aluminum

€ 8.80

M7 hemispherical nut



D00020006A

Aluminum

€ 8.80

M5 nut for stabilization hook



D00020008A

Stainless steel

€ 2.97

M5 small nut with flange



D00010003A

Stainless steel

€ 0.44

M4 nut with flange



D00010001A

Stainless steel

€ 0.44

External Fixation

UNIVERSAL SYSTEM

Carbon rod

Carbon rod
Ø 5.0 mm
Length 50 mm

B00010010A

Carbon composite

€ 2.75

Carbon rod
Ø 5.0 mm
Length 75 mm

B00010011A

Carbon composite

€ 3.03

Carbon rod
Ø 5.0 mm
Length 100 mm



B00010001A

Carbon composite

€ 3.30

Carbon rod
Ø 5.0 mm
Length 125 mm

B00010012A

Carbon composite

€ 3.85

Carbon rod
Ø 5.0 mm
Length 150 mm

B00010002A

Carbon composite

€ 4.40

Carbon rod
Ø 5.0 mm
Length 200 mm

B00010003A

Carbon composite

€ 5.50

Carbon rod
Ø 8.0 mm
Length 150 mm

B00010004A

Carbon composite

€ 5.50

Carbon rod
Ø 8.0 mm
Length 200 mm



B00010005A

Carbon composite

€ 6.60

Carbon rod
Ø 8.0 mm
Length 250 mm

B00010006A

Carbon composite

€ 7.70

External Fixation

UNIVERSAL SYSTEM

Washers

Round washer
7.5x14 mm

R00020002A

Stainless steel

€ 0.22

Round washer
5.5x15 mm



R00020003A

Stainless steel

€ 0.22

**Rectangular washer
5.5 x 15 mm**

R00030005A

Stainless steel

€ 1.10



**Rectangular washer
8.5 x 15 mm**

R00030006A

Stainless steel

€ 1.10

**Rectangular
self-centering washer
5.5x15 mm**

R00030008A

Stainless steel

€ 2.20



**Rectan gular
self-centering washer
8.2x15 mm**

R00030009A

Stainless steel

€ 2.20

**Hemispheric washer
10x16 mm**

R00030001A

Stainless steel

€ 3.30



**Hemispheric washer
8x14 mm**

R00030002A

Stainless steel

€ 3.30



**Radiolucent
hemispherical washer
10x15 mm**

R00030003A

Technopolymer

€ 4.40

**Toothed ss washer
hole 3 mm**



R00020006A

Stainless steel

€ 0.11

External Fixation

UNIVERSAL SYSTEM

Hybrid fixation replacement elements

**Pin-holder cylinder 1
hole**



C00020001A

Aluminum

€ 11.00

**Pin-holder cylinder 2
holes**



C00020002A

Aluminum

€ 15.40

**Pin-holder cylinder 3
holes**



C00020003A

Aluminum

€ 19.80

Tensioner head filler



D00030001A

Aluminum

€ 24.20

**Clamp extension for
dynamizable fixator**



E00010001A

Stainless steel

€ 12.10

**Dynamizable linear
fixator 10 cm**



F00030001B

Stainless steel

€ 140.80

**Dynamizable linear
fixator 15 cm**



F00030002B

Stainless steel

€ 150.70

**Locking hook for pins
Ø 3 mm
Length 16.5 mm x 5.0 mm
for Polilock fixator**



G00010001A

Stainless steel

€ 3.30

**Long locking hook for pins
Ø 3 mm
Length 24.5 mm x 5.0 mm
for Polilock f ixator**



G00010002A

Stainless steel

€ 3.30

**Locking hook for bars
Ø 3 mm
Length 16.5 mm x 7.5 mm
for Polilock fixator**



G00010004A

Stainless steel

€ 3.30

**Cardan joint, M5
thread**



G00020005A

Stainless steel

€ 121.00

**Transarticular round
joint**



G00020002A

Stainless steel

€ 35.20

**Linear/circular
connector**



M00010002A

Stainless steel

€ 17.60

Wire/pin clamp



M00010003A

Aluminum

€ 9.90

**Lightweight clamp for
dynamizable fixator**



M00010018A

Aluminum

€ 29.70

**Support polilock
Ø 20 mm
with 3 aligned holes**



M00010020A

Technopolymer

€ 11.00

**Support polilock
Ø 20 mm
with 3 staggered holes**



M00010021A

Technopolymer

€ 11.00

**Cap for pins
Ø 3.00 mm**



P00040002A

Rubber

€ 0.17

**Cap for pins
Ø 4.00 mm**



P00040003A

Rubber

€ 0.17

Cap for pins
Ø 1.6 mm



P00040004A

Rubber

€ 0.17

Cap for pins
Ø 1.8 mm

P00040005A

Rubber

€ 0.17

Threaded locker for pin
Length 8 mm

V00020001A

Stainless steel

€ 0.44

Threaded locker for pin
Length 12 mm



V00020002A

Stainless steel

€ 0.50

Threaded locker for pin
Length 16 mm

V00020003A

Stainless steel

€ 0.55

Hexagonal head
perforated screw
Length 21.5 mm
for K wires



V00030001A

Stainless steel

€ 5.50

Hexagonal head
perforated screw
Length 13.5 mm
for K wires

V00030002A

Stainless steel

€ 5.50

Malleable bar
Ø 3.0 mm
Length 200 mm



B00010008A

Stainless steel

€ 2.20

Screw hexagonal head
M3x20 mm



V00020013A

Stainless steel

€ 0.28

Dynamometric wire
tensioner



U00010001A

Stainless steel

€ 440.00

Hexa gon wrench
3.0 mm
with spherical tip



U00010002A

Stainless steel

€ 2.75

Wrench for hooks



U00010006A

Stainless steel

€ 44.00

Socket wrench
5.5 mm
straight handle



U00010014A

Stainless steel

€ 35.20

Double-purchase post,
M6



M00010026A

Aluminum

€ 19.80

Autoclavable container for linear radiolucent rail external fixator kit



C00040025A

Thermoplastic material

€ 33.00

Autoclavable container for dynamizable fixator kit



C00040010A

Thermoplastic material

€ 24.20

Autoclavable container for mini hybrid kit



C00040018A

Thermoplastic material

€ 33.00

Autoclavable container for linear radiolucent midi kit



C00040016A

Thermoplastic material

€ 33.00

Autoclavable container for linear aluminum midi kit



C00040013A

Thermoplastic material

€ 33.00

Autoclavable container for maxi hybrid kit



C00040015A

Thermoplastic material

€ 33.00

Autoclavable complete container for linear radiolucent mini kit



C00040008A

Thermoplastic material

€ 24.20



Porous titanium cylinders

POROUS TITANIUM CYLINDER

Porous cylinder
20 mm - Ø 8 mm



C00100001A

Titanium

€ 77.00

Porous cylinder
20 mm - Ø 12 mm



C00100002A

Titanium

€ 77.00

Porous cylinder
20 mm - Ø 16 mm



C00100003A

Titanium

€ 77.00

All the cylinders can be stabilized to a plate by means of 1.5 mm screw



Instruments

INSTRUMENTS

**Distractor standard
16 cm**



S00040010A

Stainless steel

€ 198.00

**Power distractor
17 cm**



S00040008A

Stainless steel

€ 275.00

Titan distractor



S00040011A

Aluminum

€ 1,012.00

**Guide for X-Porous TTA
osteotomy**



S00050002A

Stainless steel

€ 79.20

**Fiber-handle
screwdriver for
hexagonal screws
Ø 1.5 mm**



U00010055A

Fiber-handle

€ 90.20

**Fiber-handle
screwdriver for
hexagonal screws
Ø 2.5 mm**

U00010056A

Fiber-handle

€ 90.20

Mini Screwdriver for cross head screws interchangeable tip



U00010043A

Stainless steel

€ 202.20

Screwdriver for hex head locking screws



U00010033A

Plastic-handle

€ 60.50

Torque Limiting Screwdriver for hex head locking screws



U00010035A

Stainless steel

€ 346.50

Hexagonal quick coupling for screwdriver hexagonal screws 1.5 mm



P00070001A

Stainless steel

€ 35.20

Hexagonal quick coupling for screwdriver hexagonal screws 2.5 mm



P00070002A

Stainless steel

€ 35.20

Screwdriver with quick coupling (with P00070001A and P00070002A)



U00010030A

Fiber-handle

€ 151.80

Impactor Length 130 mm



B00010014A

Plastic

€ 28.60

Extractor hook Ø 1.0



G00050001A

Stainless steel

€ 13.20

Extractor hook Ø 1.2



G00050002A

Stainless steel

€ 13.20

Radiographic Reference Ball



S00030003A

Technopolymer

€ 113.30

**Drill sleeve Ø 2.0 mm
for locking screws,
stainless steel, with
lock**



C00010014A

Stainless steel

€ 38.50



Tutors

BODICES FOR POSTOPERATIVE MANAGEMENT OF SHOULDER AND ELBOW PROCEDURE

Bodice for front limb surgical protection 1° dx



C00090017A

85% C otton
15% Poly ester

€ 108.90

Bodice for front limb surgical protection 1° sx

C00090018A

85% C otton
15% Poly ester

€ 108.90

Bodice for front limb surgical protection 2° dx



C00090019A

85% C otton
15% Poly ester

€ 108.90

Bodice for front limb surgical protection 2° sx

C00090020A

85% C otton
15% Poly ester

€ 108.90

Bodice for front limb surgical protection 3° dx



C00090021A

85% C otton
15% Poly ester

€ 108.90

Bodice for front limb surgical protection 3° sx

C00090022A

85% C otton
15% Poly ester

€ 108.90

**Bodice for front limb
surgical protection
4° dx**

C00090023A

85% Cotton
15% Polyester

€108.90

**Bodice for front limb
surgical protection
4° sx**



C00090024A

85% Cotton
15% Polyester

€108.90

**Bodice for front limb
surgical protection
5° dx**

C00090025A

85% Cotton
15% Polyester

€121.00

**Bodice for front limb
surgical protection
5° sx**



C00090026A

85% Cotton
15% Polyester

€121.00

**Bodice for front limb
surgical protection
6° dx**

C00090027A

85% Cotton
15% Polyester

€121.00

**Bodice for front limb
surgical protection
6° sx**



C00090028A

85% Cotton
15% Polyester

€121.00

**Bodice for front limb
surgical protection
7° dx**

C00090029A

85% Cotton
15% Polyester

€121.00

**Bodice for front limb
surgical protection
7° sx**

C00090030A

85% Cotton
15% Polyester

€ 121.00

**Bodice for front limb
surgical protection
8° dx**



C00090031A

85% Cotton
15% Polyester

€ 121.00

**Bodice for front limb
surgical protection
8° sx**

C00090032A

85% Cotton
15% Polyester

€ 121.00



K wires and self-tapping threaded pins

K wires

HALF TIP K WIRES

**K wire half tip,
Ø 1.0 mm
Length 220 mm**

F00020157A

Stainless steel

€ 4.40

**K wire half tip,
Ø 1.2 mm
Length 250 mm**



F00020158A

Stainless steel

€ 4.51

**K wire half tip,
Ø 1.5 mm
Length 250 mm**

F00020159A

Stainless steel

€ 4.62

K wires

DOUBLE TROCAR K WIRES TIP

**K wires double Trocar
tip Ø 1.0 mm
Length 300 mm**



F00020178A

Stainless steel

€ 3.52

**K wires double Trocar tip
Ø 1.2 mm
Length 300 mm**

F00020179A

Stainless steel

€ 3.52

K wires double Trocar
tip Ø 1.5 mm
Length 300 mm



F00020180A

Stainless steel

€ 3.52

K wires double Trocar
tip Ø 1.8 mm
Length 300 mm



F00020181A

Stainless steel

€ 3.52

k WIRES

K WIRE NEGATIVE THRAD

K wire Ø 1.2 mm
Length 120 mm,
negative 10 mm thread



F00020182A

Stainless steel

€ 4.73

K wire Ø 1.6 mm
Length 120 mm,
negative 10 mm thread



F00020183A

Stainless steel

€ 4.73

Drill bit

DRILL BITS

Drill bit
Ø 1.5 mm
Length 100 mm

P00060005A

Stainless steel

€ 15.95

Drill bit
Ø 1.8 mm
Length 100 mm

P00060001A

Stainless steel

€ 15.95

Drill bit
Ø 2.0 mm
Length 120 mm



P00060002A

Stainless steel

€ 15.95

Drill bit
Ø 2.5 mm
Length 130 mm

P00060003A

Stainless steel

€ 15.95

Drill bit
Ø 3.2 mm
Length 130 mm

P00060004A

Stainless steel

€ 15.95

SELF-TAPPING THREADED PINS

Self-tapping threaded pins
Ø 2.0
Length 100 mm
15 mm thread

F00020173A

Stainless steel

€ 11.55

Self-tapping threaded pins
Ø 2.4
Length 100 mm
15 mm thread

F00020174A

Stainless steel

€ 11.55

Self-tapping threaded pins
Ø 2.7
Length 120 mm
15 mm thread



F00020175A

Stainless steel

€ 11.55

Self-tapping threaded pins
Ø 3.5
Length 130 mm
25 mm thread

F00020176A

Stainless steel

€ 11.55

Self-tapping threaded pins
Ø 4.5
Length 130 mm
25 mm thread

F00020177A

Stainless steel

€ 11.55

